HOMEDICS®

NO TOUCH INFRARED THERMOMETER

Le manuel in francais
commence a la page 23

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Model: NT17
P/N: TIE-210

TABLE OF CONTENTS

Introduction	3
Important Information Before Use	5
Product Identification.....	7
Description of LCD Display	8
Battery Installation	9
Switching Between Fahrenheit or Celsius ..	10
Switching between 4 Kinds of Scan Modes.....	10
Tips for Measuring Human Temperature	11
Measuring Human Temperature	12

Measuring Object/Liquid Temperature.....	13
Memory Function	14
Care and Maintenance	15
Applied Standards.....	15
Error Codes	16
Technical Specification	17
EMC Tables.....	18
Warranty.....	22

INTRODUCTION

Utilizing infrared technology, this thermometer takes temperatures in seconds by measuring heat generated by the surface skin of the forehead or other objects. This product conforms to the provisions of the EC directive MDD (93/42/EEC). Please note that this product is not intended to serve as a substitute for the advice of a physician or medical professional.

Features include:

- 6 Device Functions
Human Body/Object/Night Mode Option/
Fever Alarm/30 User Memory/Backlight
- Night Mode Option
Switch to Night mode to reduce the interference of buzzer during your baby sleeping.
- Red LED Light for Fever Alarm
This device has the LED light to remind user their readings are over 100.4°F (38°C) and plus the beep sound only in Human Body Mode.
- 30 User Memory

- Illuminated Backlight Display
- °F/°C Switch-able Function
- One-second Reading
- Auto power off for power saving
- Low-battery Indicator
Indications for battery condition and measuring range.
- Large LCD Display
- Economic design and convenience:
This is a “Non-contact” medical thermometer that enables temperature readings, designed for sanitary conditions, cleanliness, and convenience. Simply move the thermometer close to the subject's forehead or object at the distance indicated by the device.
- Instant Measurement
By using our unique technology, users can get their precise body temperature instantly and accurate.


IMPORTANT INFORMATION BEFORE USE

When using your thermometer, basic precautions should always be followed. Please read and follow all instructions and warnings before using this product. Save these instructions for future reference.


- Do not disassemble, repair, or remodel the thermometer.
- Be sure to clean the thermometer lens after each use, refer to the care and maintenance section.
- Avoid direct finger contact with the lens.
- No modification of this equipment is allowed.
- It is recommended that user may take 3 temperatures. If they are different, use the highest reading.
- Do not expose the thermometer to extreme temperature, very high humidity, or direct sunlight.
- Avoid extreme shock or dropping the device.
- Before the measurement, patients and thermometer should stay in steady state room condition for at least 30 minutes.

- Avoid measuring temperature in 30 minutes after exercise, bathing, or returning from outdoor.
- To protect the environment, dispose of empty batteries at appropriate collection sites according to national or local regulations.
- It is ill-advised to disassemble the thermometer.
- Please use the thermometer solely for its intended purpose.
- Carefully hold the device when in use to avoid dropping the device.
- Allow one minute between successive measurements as slight variations may occur if measurements are taken over a short period of time. Use average temperatures instead.
- There are no absolute body temperature standards. Keep reliable records of your personal temperature to serve as a reference for judging a fever.
- Under any circumstances, the temperature taking result is **ONLY** for reference. Before taking any medical action, please consult your physician.

PRODUCT IDENTIFICATION

DESCRIPTION OF LCD DISPLAY


BATTERY INSTALLATION

Low Battery Warning:

When the battery power becomes low, the low battery symbol “b” will appear on the display. The thermometer can still be used during this time, but the batteries should be replaced as soon as possible. If the batteries run out completely, “Lo” will be displayed along with the low battery symbol. In this case, the batteries will need to be replaced before using the thermometer again.

Replacing the Battery:

1. Gently slide the battery cover back.
2. Carefully remove the old batteries and properly discard.
3. Insert new batteries (Two 1.5V alkaline AAA Size) according to the proper polarity.
4. Slide the battery cover back on.

NOTE: Battery-operated

1. Please properly dispose of the batteries away from small children and heat.
2. It is recommended to remove the batteries if the unit will not be used for an extended period of time.
3. For long durations of non-operation, please remove all batteries from the device.
4. Batteries must be disposed of in accordance with local environmental and institutional policies.
5. Dispose of used batteries in accordance with the applicable legal regulations. Never dispose of batteries in the normal household waste.

SWITCHING BETWEEN FAHRENHEIT OR CELSIUS

Your thermometer can display results in either degrees Celsius ($^{\circ}\text{C}$) or degrees Fahrenheit ($^{\circ}\text{F}$). To switch between Celsius and Fahrenheit, while the unit is on press and hold both the Human/Object button and the Memory button for approximately 3 seconds. This will change the mode to either $^{\circ}\text{C}$ or $^{\circ}\text{F}$. Once the thermometer beeps, at which time the mode has been changed and you may release the buttons.

SWITCHING BETWEEN 4 KINDS OF SCAN MODE

1. Under power on status, you can press the Human/Object button to switch different scan mode. There are 4 kinds of mode which including Human, Object, Human/Night, and Object/Night mode (in order).
2. The beep sounds will be closed when your choice in Human/Night mode, and Object/Night mode, and the Moon symbol will appears on the LCD in both Night mode.

Human Mode

Object Mode

Human Mode
& Night Mode

Object Mode
& Night Mode

NOTE:

Each press will comes with a beep sound to ensure the setting is activated. (Except both Night modes).

TIPS FOR MEASURING HUMAN TEMPERATURE


Bear in mind that the thermometer needs to have been in the room in which the measurement is taken for at least 30 minutes before use.

- Attempting to take temperature readings from sites on the body other than the forehead may produce inaccurate results.
- The patient should remain still while the reading is being taken.
- Infrared forehead temperature readings are equivalent to oral temperature readings. In all of these cases, please consult your doctor.
- Readings taken while asleep should not be compared directly to readings taken while awake, as body temperature while asleep is typically lower.
- Do not take body temperature readings within 30 minutes of being outdoors, exercising or bathing

MEASURING HUMAN TEMPERATURE

1. Press the Power button to turn the thermometer on. The unit will run a self-test and all symbols on the display will momentarily appear.
2. Ensure that the thermometer is in Human mode; the Human symbol will be on the display. To alternate between modes press and release the Human/Object Mode button until you see the desired measurement symbol on the display.
3. Press and hold the "Scan & measure button", aiming the Measurement Sensor at the person's forehead, holding the thermometer approximately 0.5~3 cm. But don't touch the forehead.
4. Release the "Scan & measure button". In about a second, you will hear a short beep means this temperature reading has been completed and accompany with a back-light.
5. If the temperature measurement is below 38°C (100.4°F), a "Smiling Face 😊" will appear next to the reading. If the reading is 38°C (100.4°F) or above, a "Frowning Face ☹️" will be displayed and the RED LED light up.


Refer to
instruction

Fever Alarm: Only in Human Mode and Human/Night Mode

6. After about 30 seconds after use, the thermometer will automatically beep and shut off.

MEASURING OBJECT/LIQUID TEMPERATURE

1. Press the Power button to turn the thermometer on. All symbols on the display will momentarily appear.
2. Ensure that the thermometer is in Object mode; the Object symbol will be on the display. To alternate between modes press and release the Human/Object Mode button until you see the desired measurement symbol on the display.
3. Press and hold the “Scan & measure button”, aiming the Measurement Sensor at the object, holding the thermometer approximately 0.5~3 cm.
4. Release the “Scan & measure button” and the temperature reading will be displayed.
5. After about 30 seconds after use, the thermometer will automatically beep and shut off.


MEMORY FUNCTION

Memory Recall:

You can recall up to 30 measurements currently stored in memory to share with your physician or trained healthcare professional.

1. When the device is on, press once briefly on the "MEM button", then pass it again to show the last measurement accompanied by "☑ 1" symbol.
2. The "👤" symbol or "📶" symbol will appear with each measurement stored in memory to indicate whether a person or object temperature was taken.
3. Each press of the same button recalls a previous measurement, so "☑ 30" then all the way to "☑ 1".

Memory Deletion:


1. Under power on status, you may keep press MEM Button for more than 3 seconds to delete all the readings.
2. You will see "--" on display and "sound four short beeps", that mean all the memories are cleared.
3. Automatically on the 31st measurement: when the 30 memories have been used up, any new measurement will be recorded with "☑ 0" and the oldest memory deleted without you having to do anything.

All the readings will be cleared no matter record in Human mode or Object mode.

CARE AND MAINTENANCE

- **Lens/ Measurement Sensor:**
Gently clean with an alcohol swab. Do not use water to wash the thermometer lens directly.
- **Thermometer:**
Clean with a soft, dry cloth. Do not use water to rinse the device.

APPLIED STANDARDS


This product conforms to the provisions of the EC directive MDD (93/42/EEC). The following standards apply to design and/or manufacture of the products:

- **ASTM E1965-98**
Standard Specification for Infrared Thermometers for Intermittent Determination of Patient Temperature.
- **ISO 14971**
Medical devices-Application of risk management to medical devices. The Classification according to IEC/EN 60601-1 sub-clauses 5:
 - Internally powered equipment
 - IPX0
 - Equipment not suitable for use in the presence of a flammable anaesthetic mixture with air or with oxygen or nitrous oxide
 - Continuous operation

	Refer to instruction
	Type BF equipment (Sensor)
	The device shall be disposed in accordance with national laws after their useful lives
	Manufacturer's name and address

ERROR CODES

When a malfunction or incorrect temperature measurement occurs, an error message will appear as described below.

LCD DISPLAY	CAUSE	SOLUTION
	The temperature measured is higher than 1. Human thermometer mode: 43°C (109.4°F) 2. Object temperature mode: 100°C (212.0°F)	Operate the thermometer only between the specified temperature ranges. If necessary, clean the sensor tip. In the event of a repeated error message, contact your retailer or Customer Services.
	The temperature measured is higher than 1. Human thermometer mode: 34°C (93.2°F) 2. Object temperature mode: 0°C (32.0°F)	
	The operating temperature is not in the range from 16°C~40°C (60.8°F to 104°F)	Operate the thermometer only between the specified temperature ranges.

TECHNICAL SPECIFICATION

- **Measuring range :**

Human Body: 34°C~43°C (93.2°F~109.4°F)

Object: 0°C~100°C (32.0°F~212.0°F)

- **Calibration Accuracy:**

Human Body: 36°C~39°C (96.8°F~102.2°F): ± 0.2°C

else: ± 0.3°C

Object: ≤ 20°C : ± 1°C

> 20°C: ± 5%

- **Display resolution: 0.1°C**

- **Operating environment:**

16°C~40°C (60.8°F~104°F)

with relative humidity up to 95% (non condensing)

- **Storage/ Transportation environment:**

-20°C~+50°C (-4°F~+122°F)

with relative humidity up to 95% (non condensing)

- **Power supply:** 2 x 1.5V AAA size alkaline batteries
- **Weight:** approx. 80g (with batteries)
- **Dimensions:** approx. 128.5mm×48.83mm×38.85mm (L×W×H)
- **Operation Distance:** 0.5~3 cm

EMC TABLES

Guidance and manufacturer's declaration-electromagnetic emissions		
The TIE-210 Thermometer is intended for use in the electromagnetic environment specified below. The customer or the user of the TIE-210 Thermometer should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment guidance
RF emissions CISPR 11	Group 1	The TIE-210 Thermometer uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The TIE-210 Thermometer is suitable for use in all establishments, including domestic establishments and those directly connected to the public lowvoltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable	

EMC TABLES

Guidance and manufacturer's declaration-electromagnetic immunity			
The TIE-210 Thermometer is intended for use in the electromagnetic environment specified below. The customer or the user of the TIE-210 Thermometer should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

EMC TABLES

Guidance and manufacturer's declaration-electromagnetic immunity

The TIE-210 Thermometer is intended for use in the electromagnetic environment specified below. The customer or the user of the TIE-210 Thermometer should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Radiated RF IEC 61000-4-3	80 MHz to 2,5 GHz	3 V/m	<p>Portable and mobile RF communications equipment should be used no closer to any part of the TIE-210 Thermometer, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance $d = 1.2 \sqrt{P}$ $d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2,5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a, should be less than the compliance level in each frequency range^b.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol: </p>

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the TIE-210 Thermometer is used exceeds the applicable RF compliance level above, the TIE-210 Thermometer should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the TIE-210 Thermometer.

b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

EMC TABLES

Recommended separation distances between portable and mobile RF communications equipment and the TIE-210 Thermometer

The TIE-210 Thermometer is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the TIE-210 Thermometer can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the TIE-210 Thermometer as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz d=1.2 √P	80 MHz to 800 MHz d=1.2 √P	800 MHz to 2.5 GHz d=2.3 √P
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

1-YEAR LIMITED WARRANTY

**For service or repair,
do not return this unit
to the retailer. Contact HoMedics
Consumer Relations at:**

Email:
cservice@homedics.com

Phone:
1-800-466-3342
Business Hours:
8:30am-7pm EST
Monday-Friday

In USA Distributed by:
HoMedics USA, LLC
3000 N Pontiac Trail
Commerce Township, MI 48390
Printed in China

In Canada Manufactured for:
HoMedics Group Canada
A Division of HoMedics
USA, LLC
6460 Kennedy Road, Unit C.
Mississauga, ON L5T 2X4
Toll Free:
1-888-225-7378
Made in China

HoMedics sells its products with the intent that they are free of defects in manufacture and workmanship for a period of 1 years from the date of original purchase, except as noted below. HoMedics warrants that its products will be free of defects in material and workmanship under normal use and service. This warranty extends only to Consumers and does not extend to Retailers.

To obtain warranty service on your HoMedics product, contact a consumer relations representative for assistance. Please make sure to have the model number of the product available.

HoMedics does not authorize anyone, including, but not limited to, Retailers, the subsequent Consumer purchaser of the product from a Retailer or remote purchasers, to obligate HoMedics in any way beyond the terms set forth herein. This warranty does not cover damage caused by misuse or abuse; accident; the attachment of any unauthorized accessory; alteration to the product; improper installation; unauthorized repairs or modifications; improper use of electrical/power supply; loss of power; dropped product; malfunction or damage of an operating part from failure to provide manufacturer's recommended maintenance; transportation damage; theft; neglect; vandalism; or environmental conditions; loss of use during the period the product is at a repair facility or otherwise awaiting parts or repair; or any other conditions whatsoever that are beyond the control of HoMedics.

This warranty is effective only if the product is purchased and operated in the country in which the product is purchased. A product that requires modifications or adoption to enable it to operate in any other country than the country for which it was designed, manufactured, approved and/or authorized, or repair of products damaged by these modifications is not covered under this warranty.

THE WARRANTY PROVIDED HEREIN SHALL BE THE SOLE AND EXCLUSIVE WARRANTY. THERE SHALL BE NO OTHER WARRANTIES EXPRESS OR IMPLIED INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS OR ANY OTHER OBLIGATION ON THE PART OF THE COMPANY WITH RESPECT TO PRODUCTS COVERED BY THIS WARRANTY. HOMEDICS SHALL HAVE NO LIABILITY FOR ANY INCIDENTAL, CONSEQUENTIAL OR SPECIAL DAMAGES. IN NO EVENT SHALL THIS WARRANTY REQUIRE MORE THAN THE REPAIR OR REPLACEMENT OF ANY PART OR PARTS WHICH ARE FOUND TO BE DEFECTIVE WITHIN THE EFFECTIVE PERIOD OF THE WARRANTY. NO REFUNDS WILL BE GIVEN. IF REPLACEMENT PARTS FOR DEFECTIVE MATERIALS ARE NOT AVAILABLE, HOMEDICS RESERVES THE RIGHT TO MAKE PRODUCT SUBSTITUTIONS IN LIEU OF REPAIR OR REPLACEMENT.

This warranty does not extend to the purchase of opened, used, repaired, repackaged and/or resealed products, including but not limited to the sale of such products on internet auction sites and/or sales of such products by surplus or bulk resellers. Any and all warranties or guarantees shall immediately cease and terminate as to any products or parts thereof which are repaired, replaced, altered, or modified, without the prior express and written consent of HoMedics.

This warranty provides you with specific legal rights. You may have additional rights which may vary from state to state. Because of individual state regulations, some of the above limitations and exclusions may not apply to you.

For more information regarding our product line in the USA, please visit: www.homedics.com.

HoMedics® is a registered trademark of HoMedics, LLC.
©2020 HoMedics, LLC. All rights reserved.
IB-TIE210

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